

## Product datasheet for **TB433767**

### Casein Kinase 1 alpha (CSNK1A1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CSNK1A1 (NM_001271741), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC233767
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001271741</u> , <u>NP_001258670</u>
Synonyms:	CK1; CK1a; CKIa; HEL-S-77p; HLCDGP1; PRO2975
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001258670</u>
Locus ID:	1452
Cytogenetics:	5q32
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase
Protein Pathways:	Hedgehog signaling pathway, Wnt signaling pathway



[View online »](#)