

Product datasheet for **TB433777**

Osteopontin (SPP1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SPP1 (NM_001251830), transcript variant 5, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC233777
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Osteopontin (SPP1) Mouse Monoclonal Antibody [Clone ID: OTI5E4] (TA806784)
ACCN:	NM_001251830 , NP_001238759
Synonyms:	BNSP; BSPI; ETA-1; OPN
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001238759
Locus ID:	6696
Cytogenetics:	4q22.1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	ECM-receptor interaction, Focal adhesion, Toll-like receptor signaling pathway



[View online »](#)