

Product datasheet for **TB434033**

LTBR CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of LTBR (NM_001270987), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC234033 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | LTBR Mouse Monoclonal Antibody [Clone ID: OTI3G10] (TA805712) |
| ACCN: | <u>NM_001270987</u> , <u>NP_001257916</u> |
| Synonyms: | D12S370; LT-BETA-R; TNF-R-III; TNFCR; TNFR-RP; TNFR2-RP; TNFR3; TNFRSF3 |
| Storage: | Room Temperature, or 2-8°C for long term storage |
| Stability: | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP_001257916</u> |
| Locus ID: | 4055 |
| Cytogenetics: | 12p13.31 |
| Protein Families: | Druggable Genome, Transmembrane |
| Protein Pathways: | Cytokine-cytokine receptor interaction |



[View online »](#)