

## Product datasheet for **TB434181**

### SLC16A3 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SLC16A3 (NM_001206951), transcript variant 5, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234181
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001206951</u> , <u>NP_001193880</u>
Synonyms:	MCT-3; MCT-4; MCT 3; MCT3; MCT 4; MCT4
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001193880</u>
Locus ID:	9123
Cytogenetics:	17q25.3
Protein Families:	Transmembrane



[View online »](#)