

Product datasheet for **TB434248**

Nectin 3 (NECTIN3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PVRL3 (NM_001243288), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234248
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001243288 , NP_001230217
Synonyms:	CD113; CDW113; NECTIN-3; PPR3; PRR3; PVRL3; PVRR3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001230217
Locus ID:	25945
Cytogenetics:	3q13.13
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Adherens junction, Cell adhesion molecules (CAMs)



[View online »](#)