

Product datasheet for **TB434288**

Angiopoietin 1 (ANGPT1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ANGPT1 (NM_001199859), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234288
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Angiopoietin 1 (ANGPT1) Mouse Monoclonal Antibody [Clone ID: OT17F9] (TA804638)
ACCN:	<u>NM_001199859</u> , <u>NP_001186788</u>
Synonyms:	AGP1; AGPT; ANG1; HAE5
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001186788</u>
Locus ID:	284
Cytogenetics:	8q23.1
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein



[View online »](#)