

Product datasheet for **TB434488**

MAK CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of MAK (NM_001242385), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC234488 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | MAK Mouse Monoclonal Antibody [Clone ID: OT15A7] (TA808843) |
| ACCN: | NM_001242385 , NP_001229314 |
| Synonyms: | RP62 |
| Storage: | Room Temperature, or 2-8°C for long term storage |
| Stability: | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_001229314 |
| Locus ID: | 4117 |
| Cytogenetics: | 6p24.2 |
| Protein Families: | Druggable Genome, Protein Kinase |



[View online »](#)