

## Product datasheet for **TB434500**

### HDAC9 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HDAC9 (NM_001204144), transcript variant 6, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234500
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	HDAC9 Mouse Monoclonal Antibody [Clone ID: OT17G2] (TA805328)
ACCN:	<a href="#">NM_001204144</a> , <a href="#">NP_001191073</a>
Synonyms:	HD7; HD7b; HD9; HDAC; HDAC7; HDAC7B; HDAC9B; HDAC9FL; HDRP; MITR
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001191073</a>
Locus ID:	9734
Cytogenetics:	7p21.1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)