

## Product datasheet for **TB434538**

### **RUNX1T1 CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of RUNX1T1 (NM_001198626), transcript variant 7, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC234538
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>Target Detection Antibodies:</b>	RUNX1T1 Mouse Monoclonal Antibody [Clone ID: OT11H1] (TA504703)
<b>ACCN:</b>	<u><a href="#">NM_001198626</a></u> , <u><a href="#">NP_001185555</a></u>
<b>Synonyms:</b>	AML1-MTG8; AML1T1; CBFA2T1; CDR; ETO; MTG8; ZMYND2
<b>Storage:</b>	Room Temperature, or 2-8°C for long term storage
<b>Stability:</b>	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<u><a href="#">NP_001185555</a></u>
<b>Locus ID:</b>	862
<b>Cytogenetics:</b>	8q21.3
<b>Protein Families:</b>	Transcription Factors
<b>Protein Pathways:</b>	Acute myeloid leukemia, Pathways in cancer



[View online »](#)