

Product datasheet for **TB434732**

Sortilin (SORT1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SORT1 (NM_001205228), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234732
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001205228 , NP_001192157
Synonyms:	Gp95; LDLCQ6; NT3; NTR3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001192157
Locus ID:	6272
Cytogenetics:	1p21.3-p13.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Lysosome, Neurotrophin signaling pathway



[View online »](#)