

Product datasheet for **TB434775**

FBXL11 (KDM2A) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KDM2A (NM_001256405), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234775
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001256405</u> , <u>NP_001243334</u>
Synonyms:	CXXC8; FBL7; FBL11; FBXL11; JHDM1A; LILINA
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001243334</u>
Locus ID:	22992
Cytogenetics:	11q13.2
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)