

Product datasheet for **TB434846**

MAP3K13 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MAP3K13 (NM_001242317), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC234846
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MAP3K13 Mouse Monoclonal Antibody [Clone ID: OT15A2] (TA812277)
ACCN:	NM_001242317 , NP_001229246
Synonyms:	LZK; MEKK13; MLK
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001229246
Locus ID:	9175
Cytogenetics:	3q27.2
Protein Families:	Druggable Genome, Protein Kinase, Transcription Factors
Protein Pathways:	MAPK signaling pathway



[View online »](#)