

Product datasheet for **TB434988**

R Cadherin (CDH4) CytoSection

Product data:

| | |
|---------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CytoSections |
| Description: | Transient overexpression of CDH4 (NM_001252338), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC234988 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | R Cadherin (CDH4) Mouse Monoclonal Antibody [Clone ID: OTI6H5] (TA804816) |
| ACCN: | NM_001252338 , NP_001239267 |
| Synonyms: | CAD4; R-CAD; RCAD |
| Storage: | Room Temperature, or 2-8°C for long term storage |
| Stability: | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_001239267 |
| Locus ID: | 1002 |
| Cytogenetics: | 20q13.33 |
| Protein Families: | Transmembrane |
| Protein Pathways: | Cell adhesion molecules (CAMs) |



[View online »](#)