

Product datasheet for **TB435004**

NFkB p100 / p52 (NFKB2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NFKB2 (NM_001261403), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC235004
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001261403</u> , <u>NP_001248332</u>
Synonyms:	CVID10; H2TF1; LYT-10; LYT10; NF-kB2; p49/p100; p52; p100
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001248332</u>
Locus ID:	4791
Cytogenetics:	10q24.32
Protein Families:	Transcription Factors
Protein Pathways:	MAPK signaling pathway, Pathways in cancer



[View online »](#)