

Product datasheet for **TB435062**

MAN2C1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MAN2C1 (NM_001256496), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC235062
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001256496 , NP_001243425
Synonyms:	MAN6A8; MANA; MANA1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001243425
Locus ID:	4123
Cytogenetics:	15q24.2
Protein Families:	Druggable Genome
Protein Pathways:	Other glycan degradation



[View online »](#)