

## Product datasheet for **TB435304**

### PTPRZ1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PTPRZ1 (NM_001206839), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC235304
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001206839</u> , <u>NP_001193768</u>
Synonyms:	HPTPZ; HPTPzeta; phosphacan; PTP-ZETA; PTP18; PTPRZ; PTPZ; R-PTP-zeta-2; RPTPB; RPTPbeta
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001193768</u>
Locus ID:	5803
Cytogenetics:	7q31.32
Protein Families:	Druggable Genome, Phosphatase, Transmembrane
Protein Pathways:	Epithelial cell signaling in Helicobacter pylori infection



[View online »](#)