

Product datasheet for **TB435350**

Dicer (DICER1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DICER1 (NM_001271282), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC235350
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Dicer (DICER1) Mouse Monoclonal Antibody [Clone ID: OTI1G8] (TA801644)
ACCN:	<u>NM_001271282</u> , <u>NP_001258211</u>
Synonyms:	DCR1; Dicer; Dicer1e; GLOW; HERNA; K12H4.8-LIKE; MNG1; RMSE2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001258211</u>
Locus ID:	23405
Cytogenetics:	14q32.13
Protein Families:	Druggable Genome



[View online »](#)