

Product datasheet for **TB435933**

METTL11A (NTMT1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NTMT1 (NM_001286803), transcript variant 9, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC235933
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001286803</u> , <u>NP_001273732</u>
Synonyms:	AD-003; C9orf32; HOMET1A; METTL11A; NRMT; NRMT1; NTM1A
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001273732</u>
Locus ID:	28989
Cytogenetics:	9q34.11
Protein Families:	Druggable Genome



[View online »](#)