

Product datasheet for **TB436281**

RPA34 (RPA2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RPA2 (NM_001286076), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC236281
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	RPA34 (RPA2) Mouse Monoclonal Antibody [Clone ID: OT19A1] (TA500763)
ACCN:	<u>NM_001286076</u> , <u>NP_001273005</u>
Synonyms:	REPA2; RP-A p32; RP-A p34; RPA32
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001273005</u>
Locus ID:	6118
Cytogenetics:	1p35.3
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	DNA replication, Homologous recombination, Mismatch repair, Nucleotide excision repair



[View online »](#)