

Product datasheet for **TB436409**

SLAMF7 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SLAMF7 (NM_001282591), transcript variant 5, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC236409
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SLAMF7 Mouse Monoclonal Antibody [Clone ID: OTI3B3] (TA807414)
ACCN:	NM_001282591 , NP_001269520
Synonyms:	19A; CD319; CRACC; CS1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001269520
Locus ID:	57823
Cytogenetics:	1q23.3
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)