

## Product datasheet for **TB436453**

### Endothelin 3 (EDN3) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of EDN3 (NM_001302456), transcript variant 6, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC236453
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Endothelin 3 (EDN3) Mouse Monoclonal Antibody [Clone ID: OTI3E10] (TA808598)
ACCN:	<a href="#">NM_001302456</a> , <a href="#">NP_001289385</a>
Synonyms:	ET-3; ET3; HSCR4; PPET3; WS4B
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001289385</a>
Locus ID:	1908
Cytogenetics:	20q13.32
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)