

## Product datasheet for **TB436540**

### POLR3H CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of POLR3H (NM_001282885), transcript variant 6, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC236540
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	POLR3H Mouse Monoclonal Antibody [Clone ID: OTI4D8] (TA809432)
ACCN:	<a href="#">NM_001282885</a> , <a href="#">NP_001269814</a>
Synonyms:	C25; RPC8; RPC22.9
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001269814</a>
Locus ID:	171568
Cytogenetics:	22q13.2
Protein Families:	Transcription Factors
Protein Pathways:	Cytosolic DNA-sensing pathway, Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase



[View online »](#)