

Product datasheet for **TB436633**

NBL1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NBL1 (NM_001204086), transcript variant 5, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC236633
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001204086</u> , <u>NP_001191015</u>
Synonyms:	D1S1733E; DAN; DAND1; NB; NO3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001191015</u>
Locus ID:	4681
Cytogenetics:	1p36.13
Protein Families:	Secreted Protein



[View online »](#)