

Product datasheet for **TB436822**

DDB2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DDB2 (NM_001300734), transcript variant D1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC236822
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	DDB2 Mouse Monoclonal Antibody [Clone ID: OT12E12] (TA810712)
ACCN:	NM_001300734 , NP_001287663
Synonyms:	DDBB; UV-DDB2; XPE
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001287663
Locus ID:	1643
Cytogenetics:	11p11.2
Protein Families:	Druggable Genome
Protein Pathways:	Nucleotide excision repair, p53 signaling pathway, Ubiquitin mediated proteolysis



[View online »](#)