

Product datasheet for **TB437125**

NMNAT1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NMNAT1 (NM_001297778), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237125
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	NMNAT1 Mouse Monoclonal Antibody [Clone ID: OTI1F7] (TA502219)
ACCN:	NM_001297778 , NP_001284707
Synonyms:	LCA9; NMNAT; PNAT1; SHILCA
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001284707
Locus ID:	64802
Cytogenetics:	1p36.22
Protein Pathways:	Metabolic pathways, Nicotinate and nicotinamide metabolism



[View online »](#)