

Product datasheet for **TB437220**

PYCR1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PYCR1 (NM_001282279), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237220
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PYCR1 Mouse Monoclonal Antibody [Clone ID: OTI4F2] (TA811756)
ACCN:	NM_001282279 , NP_001269208
Synonyms:	ARCL2B; ARCL3B; P5C; P5CR; PIG45; PP222; PRO3; PYCR
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001269208
Locus ID:	5831
Cytogenetics:	17q25.3
Protein Pathways:	Arginine and proline metabolism, Metabolic pathways



[View online »](#)