

Product datasheet for **TB437317**

SLC22A17 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SLC22A17 (NM_001289050), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237317
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SLC22A17 Mouse Monoclonal Antibody [Clone ID: OTI2D6] (TA810195)
ACCN:	NM_001289050 , NP_001275979
Synonyms:	24p3R; BOCT; BOIT; hBOIT; NGALR; NGALR2; NGALR3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001275979
Locus ID:	51310
Cytogenetics:	14q11.2
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)