

Product datasheet for **TB437485**

LIAS CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LIAS (NM_001278590), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237485
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001278590</u> , <u>NP_001265519</u>
Synonyms:	HGCLAS; HUSSY-01; LAS; LIP1; LS; PDHLD
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001265519</u>
Locus ID:	11019
Cytogenetics:	4p14
Protein Pathways:	Lipoic acid metabolism, Metabolic pathways



[View online »](#)