

Product datasheet for **TB437495**

TAL1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TAL1 (NM_001290403), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237495
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TAL1 Mouse Monoclonal Antibody [Clone ID: OTI5H1] (TA503922)
ACCN:	<u>NM_001290403</u> , <u>NP_001277332</u>
Synonyms:	bHLHa17; SCL; tal-1; TCL5
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001277332</u>
Locus ID:	6886
Cytogenetics:	1p33
Protein Families:	Adult stem cells, Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors



[View online »](#)