

Product datasheet for **TB437525**

MSF (SEPT9) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of 44083 (NM_001293698), transcript variant 11, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237525
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MSF (SEPT9) Mouse Monoclonal Antibody [Clone ID: OT113B1] (TA812821)
ACCN:	NM_001293698 , NP_001280627
Synonyms:	AF17q25; MSF; MSF1; NAPB; PNUTL4; SEPT9; SeptD1; SINT1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001280627
Locus ID:	10801
Cytogenetics:	17q25.3
Protein Families:	Druggable Genome



[View online »](#)