

Product datasheet for **TB437669**

Prostatic Acid Phosphatase (ACPP) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ACPP (NM_001292037), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237669
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001292037 , NP_001278966
Synonyms:	5'-NT; ACP-3; ACP3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001278966
Locus ID:	55
Cytogenetics:	3q22.1
Protein Families:	Druggable Genome, Phosphatase, Transmembrane
Protein Pathways:	Riboflavin metabolism



[View online »](#)