

Product datasheet for **TB437749**

GDF 9 (GDF9) CytoSection

Product data:

| | |
|---------------------------------------|--|
| Product Type: | CytoSections |
| Description: | Transient overexpression of GDF9 (NM_001288827), transcript variant 5, in HEK293T cells, paraffin embedded controls for ICC/IHC staining |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC237749 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| ACCN: | <u>NM_001288827</u> , <u>NP_001275756</u> |
| Synonyms: | POF14 |
| Storage: | Room Temperature, or 2-8°C for long term storage |
| Stability: | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP_001275756</u> |
| Locus ID: | 2661 |
| Cytogenetics: | 5q31.1 |
| Protein Families: | Adult stem cells, Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Secreted Protein, Stem cell - Pluripotency, Stem cell relevant signaling - TGFb/BMP signaling pathway, Transmembrane |



[View online »](#)