

Product datasheet for **TB437896**

Carboxypeptidase B2 (CPB2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CPB2 (NM_001278541), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC237896
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Carboxypeptidase B2 (CPB2) Mouse Monoclonal Antibody [Clone ID: OT13D5] (TA807474)
ACCN:	NM_001278541 , NP_001265470
Synonyms:	CPU; PCPB; TAFI
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001265470
Locus ID:	1361
Cytogenetics:	13q14.13
Protein Families:	Druggable Genome, Protease, Secreted Protein
Protein Pathways:	Complement and coagulation cascades



[View online »](#)