

Product datasheet for **TB438231**

Cytochrome P450 26B (CYP26B1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CYP26B1 (NM_001277742), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238231
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Cytochrome P450 26B (CYP26B1) Mouse Monoclonal Antibody [Clone ID: OTI4A8] (TA811497)
ACCN:	NM_001277742 , NP_001264671
Synonyms:	CYP26A2; P450RAI-2; P450RAI2; RHFCA
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001264671
Locus ID:	56603
Cytogenetics:	2p13.2
Protein Families:	Druggable Genome, P450
Protein Pathways:	Retinol metabolism



[View online »](#)