

Product datasheet for **TB438261**

IL20 Receptor alpha (IL20RA) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of IL20RA (NM_001278723), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238261
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	IL20 Receptor alpha (IL20RA) Mouse Monoclonal Antibody [Clone ID: OTI7H7] (TA811623)
ACCN:	NM_001278723 , NP_001265652
Synonyms:	CRF2-8; IL-20R-alpha; IL-20R1; IL-20RA
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001265652
Locus ID:	53832
Cytogenetics:	6q23.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway



[View online »](#)