

Product datasheet for **TB438425**

SynCAM (CADM1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CADM1 (NM_001301043), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238425
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001301043</u> , <u>NP_001287972</u>
Synonyms:	BL2; IGSF4; IGSF4A; Necl-2; NECL2; RA175; sgiGSF; ST17; sTSLC-1; SYNCAM; synCAM1; TSLC1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001287972</u>
Locus ID:	23705
Cytogenetics:	11q23.3
Protein Families:	Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs)



[View online »](#)