

## Product datasheet for **TB438446**

### CD105 (ENG) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ENG (NM_001278138), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238446
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD105 (ENG) Mouse Monoclonal Antibody [Clone ID: OTI8A1] (TA506170)
ACCN:	<a href="#">NM_001278138</a> , <a href="#">NP_001265067</a>
Synonyms:	END; HHT1; ORW1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001265067</a>
Locus ID:	2022
Cytogenetics:	9q34.11
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane



[View online »](#)