

## Product datasheet for **TB438457**

### PP5 (PPP5C) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PPP5C (NM_001204284), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238457
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PP5 (PPP5C) Mouse Monoclonal Antibody [Clone ID: OTI2E12] (TA500585)
ACCN:	<a href="#">NM_001204284</a> , <a href="#">NP_001191213</a>
Synonyms:	PP5; PPP5; PPT
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001191213</a>
Locus ID:	5536
Cytogenetics:	19q13.32
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	MAPK signaling pathway



[View online »](#)