

Product datasheet for **TB438506**

RUNX2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RUNX2 (NM_001278478), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238506
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	RUNX2 Mouse Monoclonal Antibody [Clone ID: OT13E12] (TA813668)
ACCN:	NM_001278478 , NP_001265407
Synonyms:	AML3; CBF-alpha-1; CBFA1; CCD; CCD1; CLCD; OSF-2; OSF2; PEA2aA; PEBP2aA
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001265407
Locus ID:	860
Cytogenetics:	6p21.1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)