

Product datasheet for **TB438653**

Activin Receptor Type IIA (ACVR2A) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ACVR2A (NM_001278579), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238653
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Activin Receptor Type IIA (ACVR2A) Mouse Monoclonal Antibody [Clone ID: OT11F2] (TA807404)
ACCN:	NM_001278579 , NP_001265508
Synonyms:	ACTRII; ACVR2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001265508
Locus ID:	92
Cytogenetics:	2q22.3-q23.1
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, TGF-beta signaling pathway



[View online »](#)