

Product datasheet for **TB438840**

CIITA CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CIITA (NM_001286403), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238840
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CIITA Mouse Monoclonal Antibody [Clone ID: OTI7B12] (TA812200)
ACCN:	NM_001286403 , NP_001273332
Synonyms:	C2TA; CIITAIV; MHC2TA; NLRA
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001273332
Locus ID:	4261
Cytogenetics:	16p13.13
Protein Pathways:	Antigen processing and presentation, Primary immunodeficiency



[View online »](#)