

## Product datasheet for **TB438868**

### C16orf9 (FAM234A) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FAM234A (NM_001284497), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238868
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	C16orf9 (FAM234A) Mouse Monoclonal Antibody [Clone ID: OTI3B3] (TA505567)
ACCN:	<a href="#">NM_001284497</a> , <a href="#">NP_001271426</a>
Synonyms:	C16orf9; gs19; ITFG3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001271426</a>
Locus ID:	83986
Cytogenetics:	16p13.3
Protein Families:	Transmembrane



[View online »](#)