

## Product datasheet for **TB438960**

### Artemis (DCLRE1C) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DCLRE1C (NM_001289078), transcript variant g, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC238960
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_001289078</a> , <a href="#">NP_001276007</a>
Synonyms:	A-SCID; DCLREC1C; RS-SCID; SCIDA; SNM1C
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001276007</a>
Locus ID:	64421
Cytogenetics:	10p13
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle, Non-homologous end-joining, Primary immunodeficiency



[View online »](#)