

Product datasheet for **TB439477**

Periostin (POSTN) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of POSTN (NM_001286666), transcript variant 6, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC239477
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Periostin (POSTN) Mouse Monoclonal Antibody [Clone ID: OT12B2] (TA804575)
ACCN:	NM_001286666 , NP_001273595
Synonyms:	OSF-2; OSF2; PDLPOSTN; PN
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001273595
Locus ID:	10631
Cytogenetics:	13q13.3
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)