

## Product datasheet for **TB439613**

### Miz1 (ZBTB17) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ZBTB17 (NM_001287603), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC239613
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Miz1 (ZBTB17) Mouse Monoclonal Antibody [Clone ID: OTI7D8] (TA811343)
ACCN:	<a href="#">NM_001287603</a> , <a href="#">NP_001274532</a>
Synonyms:	MIZ-1; pHZ-67; ZNF60; ZNF151
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001274532</a>
Locus ID:	7709
Cytogenetics:	1p36.13
Protein Families:	Transcription Factors
Protein Pathways:	Cell cycle



[View online »](#)