

Product datasheet for **TB439924**

MSH6 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MSH6 (NM_001281493), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC239924
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MSH6 Mouse Monoclonal Antibody [Clone ID: OTI5D1] (TA807929)
ACCN:	NM_001281493 , NP_001268422
Synonyms:	GTBP; GTMBP; HNPCC5; HSAP; MMRC3; p160
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001268422
Locus ID:	2956
Cytogenetics:	2p16.3
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Colorectal cancer, Mismatch repair, Pathways in cancer



[View online »](#)