

## Product datasheet for **TB439925**

### MSH6 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MSH6 (NM_001281494), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC239925
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MSH6 Mouse Monoclonal Antibody [Clone ID: OT15D1] (TA807929)
ACCN:	<a href="#">NM_001281494</a> , <a href="#">NP_001268423</a>
Synonyms:	GTBP; GTMBP; HNPCC5; HSAP; MMRC53; p160
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001268423</a>
Locus ID:	2956
Cytogenetics:	2p16.3
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Colorectal cancer, Mismatch repair, Pathways in cancer



[View online »](#)