

## Product datasheet for **TB440126**

### IGF1 Receptor (IGF1R) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of IGF1R (NM_001291858), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC240126
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	IGF1 Receptor (IGF1R) Mouse Monoclonal Antibody [Clone ID: OTI4C4] (TA807835)
ACCN:	<a href="#">NM_001291858</a> , <a href="#">NP_001278787</a>
Synonyms:	CD221; IGFIR; IGFR; JTK13
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001278787</a>
Locus ID:	3480
Cytogenetics:	15q26.3
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	Adherens junction, Colorectal cancer, Endocytosis, Focal adhesion, Glioma, Long-term depression, Melanoma, Oocyte meiosis, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer



[View online »](#)