

## Product datasheet for **TB492555**

### SARS-CoV-2 N Protein CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SARS-CoV-2 N (YP_009724397) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Virus
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone VC102555
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	SARS-CoV-2 N Protein Mouse Monoclonal Antibody [Clone ID: OTI8G10] (TA814432)
ACCN:	<u><a href="#">YP_009724397</a></u>
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Locus ID:	43740575



[View online »](#)