


## Product datasheet for **TC400001**

### HEK293T Cells Un-transfected Negative Control

#### Product data:

Product Type:	CytoSections
Description:	Untransfected HEK293T cells, negative FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Control Type:	Negative
Species:	Human
Storage:	RT
Stability:	Stable for at least 1 year from receipt of products under proper storage and handling conditions.

#### Product images:



Immunohistochemistry staining on Negative Control cytosection TC400001 with rabbit anti DDK clone OTIR5G2 C/N [TA592569] at 1:2000 dilution 20m RT. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D11-18]). Positive stain shown with the brown chromogen present. 5 micron sections, 40x magnification.



[View online »](#)